


**INSTRUCTIONAL MANUAL  
CAT. 64830**

**PAMCELL**




## Introduction & Features

PAMCELL™ is a 3-dimensional cell culture plate which enables the formation of uniform and a wide range of sized spheroids. Plates are composed of hexagonally-arrayed spherical particles coated with cells adhesively controlled functional groups on the surface. These plates allow for the control of cell-to-substratum and cell-to-cell interactions. Plates are available in two types: R Series and T Series.

### Features

- Enables formation of uniform and wide range sized spheroids.
- Optically transparent for in situ microscopic observation in automated high throughput screening systems.
- Provides various patterns on a single plate without physical barriers.

## Product Information



- **Sterilization:** Ethylene oxide gas sterilization
- **Packaging material:** 2-layer packaging with Nitrogen gas sealed  
[Inner- Gas permeable bag / Outer - Aluminum coated bag]

**PAMCELL™ plates are only for research use.**

## Instructions

1. PAMCELL can be used without any pretreatment, like as a common culture plate. There is no dedicated observation equipment. It is suggested to open the package in the clean bench or under sterilized environment.
2. For 96-well plate, required media is 0.1~0.3 ml for each well (recommendation: 0.3 ml), and for 6-well plate is 4~5 ml.
3. Recommended cell seeding density (Tumor cell):
  - 1) 10,000~40,000 cells for each well of 96-well plate
  - 2) 3~12 x 10<sup>5</sup> cells for each well of 6-well plate

\*The volume of media and cell density may need to be optimized for specific cell type
4. Typically ideal period for exchanging (or partial exchanging) media is 2-3 days, although it can vary for the cell type.
5. Within 24 hours, cells will be attached on the cell plate. It is recommend not to apply strong vibration or shock on the plate during initial 24 hr. Between 48~72 hr, cells spontaneously migrate to the micropads of the plate.
6. Pipetting is the best method for harvesting spheroids from the plate. It is possible to use commercially available cell scraper for harvesting, although it is not recommended.

## Storage and Shelf Life

Keep the packages sealed and stored at room temperature (20~25°C) and out of direct sunlight. We recommend to use the PAMCELL plates at least 12 months from date of production. The production date is printed on the package. Once opened, use rapidly.

## Trouble Shooting

Problems	Possible Solutions
Relatively large difference in cell density at each micropad during initial cell adsorption period	Use a cell strainer for cell seeding
Small number of cells on a micropad	Increase cell seeding density.
Fluorescence background from the bottom plate	Use PAMCELL plate made with fluorescence free COP film as substrate.

### Reference Data

- ◇ Three 3D tumor cell spheroids cultured on PAMCELL
  - Schematic drawings are based on spheroids cultured on (100 / 200 μm)\* plate

#### FaDu (Epidermoid)

- Medium: DMEM/FBS, Cell density: 10,000 cells/well  
 - Spheroid shape: Bell type

$104 \pm 8 \mu\text{m}$   
 $84 \pm 7 \mu\text{m}$

T001\_A2 : 100 / 200 μm 2 day

200 μm

T001\_C2 : 100 / 300 μm 2 day

200 μm

T001\_C4 : 200 / 400 μm 2 day

200 μm

---

#### A549 (Lung carcinoma)

- Medium: DMEM/FBS, Cell density: 20,000 cells/well  
 - Spheroid shape: Bell type

$90 \pm 6 \mu\text{m}$   
 $40 \pm 3 \mu\text{m}$

T001\_A2 : 100 / 200 μm 2 day

200 μm

T001\_A4 : 200 / 300 μm 7 day


200 μm

T002\_C5 : 400 / 1,200 μm 7 day

200 μm

\*100 μm micropad diameter, 200 μm center-center spacing between micropads: indicates as (100/200 μm)

FaDu Epidermoid                      HT-29 (Adenocarcinoma)                      A4549 (Lung Carcinoma)




[25° tilted FE-SEM images of three tumor cells]

Cell Type	Height (H, µm)			Diameter (D, µm)			H/D ratio		
	Aver.	Min.	Max.	Aver.	Min.	Max.	Aver.	Min.	Max.
FaDu	84	58	126	104	83	154	0.79	0.69	0.82
HT-29	101	86	110	126	97	150	0.81	0.73	0.89
A549	40	29	75	90	77	126	0.43	0.38	0.60

[Dimension analysis of three tumors cells from confocal microscopic images]

- ◇ Fluorescence microscopic images: Actin and Nuclei staining of A549 cells
  - Plate: COP (Fluorescence free) / 96-well plate (250/350 µm)
  - Microscopy: Olympus IX51 – 4x/0.13, 40x/0.55 Lens



- ◇ High Content screening imaging : Calcein AM staining of A549 cells (Treated with Niclosamide)
- Plate: COP/ 96-well plate (R100)
  - Instrument: Molecular Device ImageXpress Micro4

