

PROVIDING SOLUTIONS:

TEM Embedding Resins

Embedding adds desired strength and support to tissue for thin sectioning. It also adds size to specimens for easy handling. Knowing the different properties of several embedding mediums can aid in determining what type of resin would be most suitable for the type of specimen being examined.

Desirable properties for all embedding mediums:

1. Easily available
2. Uniformity from one batch to another
3. Solubility in dehydrating agents
4. Low viscosity as a monomer
5. Uniform polymerization
6. Little volume change on polymerization
7. Ease of sectioning
8. Stability under electron beam
9. Doesn't extract cell constituents
10. Good chromatic (OLM) stainability
11. Low e-scattering = Contrast
12. Good preservation of fine structure
13. Low background granularity

The Three Main Types of Embedding Resins

1. Epoxy Resins

Epoxy resins are most widely used since they have most of the desired properties.

Epoxy : Anhydride ratio with higher ratio creating more linear short chain polymers for better cutting.

Weight per Epoxide (WPE): # grams of resin per 1 gram equivalent of epoxide.



Physical Appearance: Transparent, yellowish

Desirable Characteristics

1. Range from viscous liquids to fusible solids depending on the molecular weight
2. Low shrinkage (<2%)
3. Even polymerization
4. Relatively stable with respect to light, heat, and oxygen
5. Relatively stable under electron bombardment
6. Lose only ~25% of the mass of a section under normal operating conditions
7. Can be minimized by exposing the section first to a low intensity beam and then to a gradually more intense beam
8. Maintains three-dimensional structure of the tissue specimens.

Undesirable Characteristics

1. High viscosity
2. Causes severe irritation on prolonged or repeated contact
3. Some (small) loss of contrast between tissue and background



Alpha-amylase on Lowicryl HM20 section of rat pancreas.

Some Specific Resins

	Cat. No.
Spurr's	14300 (kit)
Low viscosity (60 cps), toxic; accelerator BDMA. ETOH soluble, poor staining, more shrinkage	
EMbed 812	14120 (kit)
Uniformity viscosity (150–200 cps) Hardness can be modified by altering DDSA : MNA ratio with a higher portion of MNA making a harder block with a 1:1 ratio providing a very satisfactory hardness. Hardeners: MNA/NMA, Flexibilizer: DDSA, Optional Plasticizer: DBP, Accelerators: DMP-30, DMP-10, BDMA	
Araldite	13920
High viscosity (502: 3,000cps, 6005: 1,400 cps) ETOH soluble	
EMbed/Araldite	13940 (kit)
An excellent sectioning and staining with lower viscosity than 100% Araldite	
Durcupan	14020 (kit 260 g)
Water soluble, low viscosity polyepoxide with very little shrinkage.	
Durcupan ACM	14040 (kit 1200 ml)
2 part polyepoxide with low viscosity and very little shrinkage	
Epo-Fix	1232 (kit)
Originally used for metallographic/materials cold mounting but very good sectioning properties, low viscosity, and very little shrinkage has made it a valuable medium for fibers, papers, powders, and other materials designated for thick and/or thin sectioning. Low exothermic RT polymerization (8 hr.) or 2 hr. at 60° C.	

2. Methacrylate (Acrylic Resins)

Methacrylate (acrylic resins) are used extensively for immunocytochemistry. They polymerize unevenly and are unstable under electron bombardment.

Physical appearance: Colorless, transparent

Desirable Characteristics

1. Large size sections easily cut
2. Stain readily with excellent contrast obtainable
3. UV polymerization — no OsO₄

Undesirable Characteristics

1. Uneven polymerization
2. Too much shrinkage (15-20%)
3. Reacts with O₂
4. Resins are powerful lipid solvents. Cannot use acetone
5. Hardened resin lacks stability under electron bombardment
6. 50% of mass (may be) lost on electron irradiation followed by flow of remaining resin resulting in distortion of macromolecular structure of tissue.
7. More difficult to section

Some Specific Resins

	Cat. No.
L.R. Gold (OLM histochemistry)	14370
LR White (EM)	
Very low viscosity, available in different hardness grades:	
Soft	14384
Medium	14381
Hard	14383
Lowicryl (EM)	
High cross-linking, low viscosity. Excellent for immunocytochemical work K4M (-35° C)	
1-step kit	14335
HM20 (-70° C)	14340
HM20 1-step kit	14345
KM11 (-60° C)	14350
HM23 (-80° C)	14360

3. Polyester Resins

Polyester resins are similar to epoxy resins. They have excellent properties but are usually less easily available, however some components are not stable during storage.

NOTE: Initiator medium contains Benzoyl peroxide and Cobalt Naphthenate (CN) accelerators which when mixed together may cause explosion.

Polymerize by: Light, heat, O₂

Physical appearance: Colorless, transparent

Technical Tip: Embedding Resin Stain Permeability


Polyester embedding resins are more permeable to acid stains. Epoxy embedding resins are more permeable to basic stains.

Reference: Bryant V. and Watson, J.H.L. (1967). *A comparison of light microscopy staining methods applied to a polyester and three epoxy resins.* Henry Ford Hosp. Med. Bull., 15:65

Some Specific Resins

	Cat. No.
Vestopal W — Stains easily	
Rigolac	
Styrene monomer	14650
Selectron — First embedment	
Beetle — Rapid embedment	

Glutaraldehyde-fixed kidney specimen



References

Glauert, Audrey M. 1991. *Fixation, Dehydration and Embedding of Biological Specimens*, pp. 123-125. Amsterdam, New York, Oxford: North-Holland Publishing Company.

Hayat, M.A. 1981. *Principles and Techniques of Electron Microscopy*, vol. 1, 2nd ed., pp. 154 -170. Baltimore: University Park Press.

Bozzola, John J. 1999. *Electron Microscopy: Principles and Techniques for Biologists*, 2nd ed, pp. 36-43