# **Challenges with Resins**

Every step in preparation of samples for TEM is the most critical step. Ensure that you give the same level of attention to each of them:

#### **Fixation Dehvdration** Infiltration **Embedment**

A sample may be well fixed and present beautiful morphological features, but if blocks are unsectionable - either too hard or not completely polymerized - or have holes in the tissue, then the sample is ruined. Here we address the most common issues surrounding epoxy resins:

- 1. Blocks that are too hard
- 2. Holes in the section
- 3. Incomplete polymerization/gummy resin
- **4.** Improper polymerization

#### **Epoxy Resin Issues**

**DMAE** 

#### 1. Blocks that are too hard

0.8%

Increase the percentage of \*Flexibilizer\* in the resin recipe.

Embed 812 Recipe		<b>DER 332-732 Recipe</b>	
812	50%	<b>DER 332</b>	45%
*DDSA*	30%	*DER 732*	20%
NMA	20%	DDSA	35%
DMP-30	1.5 – 2.0%	DMP-30	2%
Spurr's Recipe		Araldite Recipe	
Spurr's Reci	ipe	Araidite Reci	pe
Spurr's Reci ERL 4221	<b>ipe</b> 25%	*DDSA*	<b>pe</b> 60%
	•		

NOTE: Araldite 6005 is very viscous before polymerization and a very hard block after. This resin should only be used with hard samples when infiltration is not an issue.

**NOTE:** There are 3 different accelerators used in epoxy embedding resins, BDMA, DMAE, and DMP-30. Each one can be substituted for the other with BDMA having the lowest viscosity. The amount added to the final total mass are as follows:


BDMA	2.5 – 3.0%
DMP-30	1.5 – 2.0%
MAE	0.6 - 0.8%

#### 2. Holes in Sections

Holes in sections are caused by:

- Poor infiltration caused by residual water in tissue after incomplete dehydration.
- Incomplete polymerization caused by moisture in tissue, dehydrant, and/or resin.
- Large defects can be caused by residual trimming artifacts which may cause some pulling out of some cellular components.





Transmission Electron micrograph of Zika Virus. Virus particles are 40 nm in diameter, with an outer envelope and an inner dense core.

Courtesy of Cynthia Goldsmith, CD.

## **TECHNICAL TIP**

# The Storing of Embedding Resin **Mixtures**

Our experience over the years shows that the best way to store remaining embedding resin mixtures is in disposable plastic syringes.

#### **Procedure:**

Take a fresh "all plastic" syringe (EMS Cat. #72520-72529). Choose the appropriate size syringe depending on the volume mixture to be stored.

Remove the tip cover.

Slowly draw in the mixture.

Remove any air space by holding the syringe with the tip in the upward position pulling the plunger slowly to clear the resin from the tip as well as letting the air move towards the tip.

Slowly push the plunger upward until the resin appears at the end of the tip.

Put on the tip cover.

Wrap the syringe with aluminum foil, label it (name, mixture, and date) and store it in the refrigerator.

This method protects the mixture from contamination by condensation, and it is also more convenient for future use.















# **Challenges with Resins**

#### **Epoxy Resin Issues (continued)**

# 3. Incomplete Polymerization/Gummy Resin

Incomplete polymerization or gummy resin is often caused by the presence of residual dehydrant in the tissue, especially when ETOH is used. To avoid this, utilize propyene oxide after ETOH dehydrant for resin infiltration. Glass distilled acetone is often used but can act as a scavenger. Finally, ensure complete infiltration with 100% resin.

For successful embedment follow the following hints:


- Be sure to use fresh or securely sealed bottles of dehydrants.
- Use the purest dehydrants for the final dehydration and infiltration steps: 200 proof ETOH (not denatured) and glass distilled acetone or propylene oxide.
- · Be aware of relative humidity in your lab. All dehydrants and resins are hygroscopic and will quickly absorb moisture.
- · Be very accurate when measuring the accelerator for your resin, as excess accelerator will make the block harder and brittle.

### Methacrylates/Acrylic Resin Issues

# Improper Polymerization

For Methacrylates or Acrylic Resins, such as LR White, Gold, or Lowicryl

- Acrylic resins react with O2 so must be polymerized using gelatin capsules or PTFE molds covered with parafilm or some material which is non-permeable to gases. NOTE: DO NOT use BEEM capsules, as they are gas-permeable!
- Acetone can not be used for dehydration when acrylic resins are to be used.
- · Polymerization is an exothermic reaction and left exposed will generate enough heat to compromise immunocytochemical reactions. Polymerize in 50-60°C oven or UV.



Transmission electron microscope image of a region in the Drosophila germarium. Specific cells in the germarium contain synaptonemal complex in their nuclei.

EMS Catalog supplies mentioned	Cat. No.
Kits	
EMbed 812	14120
Araldite EMbed 812	13940
Spurr's	14300
DER 332-732	14000
Araldite 6005	13920
Lowicryl K4M	14330
LR White (Medium Grade)	14380
All Plastic Syringes	72520-72529

Transmission electron microscope image of synaptonemal complex in nuclei of Drosophila germarium cell.

Louisa Howard, Dartmouth College.

