Endless Possibilities ...

This Kirsch Note (KN) provides a general outline of steps taken to chemically prepare samples for TEM and SEM. If there is a number in parentheses at the end of a step, that is a hyperlink to the KN that contains more information on that particular step. Links to the KNs are also listed at the end of the article.

rsch

Protocols:

Start to Finish

Bio EM Processing –

Bio EM Processing

ectron Mic

© 2018

Chemical Processing – TEM Embedment (#2A) – SEM Drying (#2B)

NOTE: Follow all general lab safety protocols. (#1)

Microscopy Academy

- 1. Samples for TEM, preferably perfused, are placed in a few drops of fixative and cut into small cubes 1 mm³. If the samples have an orientation, such as skeletal muscle, they should be cut with an aspect ratio such that the narrower edge will form the block face. (#3) (#4)
- 2. Cut samples should be placed in a suitable vial containing the fixative of choice.
- 3. The samples are then washed three times in the same 0.1M buffer as used in the fixative.
- 4. Next the samples are immersed in OsO^4 , 2-4% in aqueous or buffer.
- 5. The samples are wash three more times in DI water.
- 6. Samples are taken through a graded series of dehydrant, ETOH or acetone, typically starting at 50%, 70%, 80%, 95%, and finally three changes of 100%. **NOTE:** If using ethanol, do not use denatured alcohol. (#5)
- NOTE: For SEM, stop here and refer to HMDS (#2B) OR CPD (#6A) (#6B)
- 7. The samples are now infiltrated with the final embedment resin of choice (see #7A) starting with a 2:1 ratio of final dehydrant to resin followed by a change to 1:1 ratio and finally 2 changes of 100% resin. (#7A) (#7B)
- 8. Final embedment is done in the same 100% resin as infiltration.
- 9. Desired embedment mold is selected, a drop or two of resin is put in, the sample placed near the tip, additional resin added to fill, and label is inserted. (#8)
- 10. Embedment trays are placed in 70° C oven overnight for polymerization.
- 11. Polymerized blocks are removed for the oven and brought to RT to be removed from mold.

Sectioning – Thick and Thin

- 1. The blocks are removed from the mold and placed in the ultramicrotome chuck for trimming.
- 2. A razor blade is used to remove excess resin from the top of the block face until tissue is encountered.
- 3. The sample should then be faced off with an old diamond or a glass knife to provide a flat reflective surface. (#9)
- 4. The ROI is identified and excess materials is removed to obtain a rectangle or trapezoid shaped block face $\leq 2 \text{ mm}^2$.
 - 5. The block/chuck assembly is placed in the ultramicrotome arm with the specimen arc adjustment in the vertical position and the degree numbers to the right.
 - A knife is placed in the knife holder and clearance angle set 4° for glass and 6° for Diatome diamond knives.
 - 7. The knife assembly is advanced manually to a position close to the sample and locked.
 - 8. Observing the block face/knife position the block is rotated to align the bottom edge of the block parallel to the knife edge.

- 9. The block arc and knife arc are adjusted so a complete section, top to bottom and left to right respectively, will be cut.
- Using a knife with a boat, add water to the boat and cut sections $0.5 2 \mu m$ thick. 10.
- 18 Electron Sections are retrieved from the water, placed on a glass slide, excess water is evaporated and section is stained. (#10A RMC) (#10B Leica)
 - The ROI is selected, and the block is trimmed accordingly if needed. (#9)
 - The block/knife assembly is placed back in the microtome and a thin section knife is 13. inserted into the knife holder, secured, and clearance angle is set.
 - 14. The same process as steps 7-10 are performed, except section thickness is 60-100 nm. (#11)
 - 15. Thin sections are put onto grids for post staining. (#12)

Post Staining

4

- Grids are selected to be stained. (#13) 1.
- 2. A piece of parafilm is taped to a flat surface and a Petri dish cover to cover the film is selected.
- 3. A drop of fresh UA for each grid is placed on the parafilm and a grid, with section side down, is floated on top of the drop.
 - The Petri dish cover is placed over the drops with grids and a dark cover is placed over the entire area for the time desired for staining.
- 5. Three beakers are filled with DI water and each grid is dipped 10 or more times in each beaker.
- 6. Grids are dried and are now ready for lead staining.
- 7 Another piece of parafilm is taped to the surface and 3-4 pellets of NaOH are placed on the film and covered with Petri dish cover.
 - Fresh filtered drops of lead stain are placed on the film and pre-stained grids are floated section side down on the lead drops for desired duration.
- ectron Micro8. Grids are removed and dipped 10 or more times in fresh DI water.
 - 10. Grids are dried and ready to observe in the TEM. (#13)

TEM Observation (#14) This KN covers general problem solving for TEM observation.

Linked List of KNs Mentioned in **Body of Text** © 2018 Elec

lectron Microscol

- General Laboratory Safety #1.
- #2A. **Biological Processing Steps for TEM**
- 2B. SEM Biological HMDS Processing
- #3. Fixatives for Electron Microscopy
- #4. Buffers for Electron Microscopy
- #5. Dehydrants for Electron Microscopy
- #6A. Critical Point Drying with the EMS 3100
- 6B. Critical Point Drying with EMS Q850 Electron
- #7A. **TEM Embedding Resins**
- 7B. Challenges with Resins
- #8. Embedding Samples
- Block Trimming for TEM **#9**...ce
- #10A. Thick Sectioning for TEM using the RMC Power Tome PC

2

8+ (0)

- 10B. Thick Sectioning using the Leica UC7 Ultramicrotome
- #11. How to Thin Section TEM Samples
- #12. Three Techniques for Retrieving Thin Sections
- Preparation of Post Staining Solutions #13.0
- #14. **TEM Imaging Issues and Solutions**

Electron Microscopy Sciences | 1560 Industry Road | Hatfield, PA 19440 P 215-412-8400 | F 215-412-8450 | info@emsdiasum.com