PROTOCOL:

Biological Processing Steps for TEM

Chemical processing of biological tissues for microscopic examination has evolved to keep pace with increased levels of detail seen with newer technologies. From simple formalin fixation, paraffin embedment and 0.5 mm sections for OLM into a multichemical, epoxy resin embedment and 60 nm sections for TEM. Even this well-established TEM process has evolved further with the use of microwaves to significantly decrease the time and potential processing artifacts. This generalized procedure gives the acceptable changes, concentrations and times for both bench top and microwave processing.

Setup requirements


- Program microwave for desired processing times.
- Get ORGANIZED!! Have equipment and solutions ready.
- Samples should be placed on a rotator when using Benchtop processing.
- Setup vacuum chamber and/or agitation, if needed.

All of the following steps can be carried out in the 1.7 ml microfuge tubes, scintillation vials, or Petri dishes.

NOTE: Acetone or propylene oxide (PO) cannot be used for dehydration if plastic Petri dishes are used.

Step			Temp.	Microwave	Lynx II/Benchtop
1.	Initial fixation (Karnovsky's)		37° C	2:30 min	2 hr.
2.	Buffer rinse 3 changes		37° C	60 sec. ea.	10 min. ea.
3.	2-4% OsO4 in DI water		37° C	2:30 min	2 hr.
	(Sometimes 2% Pot	ganate in DI is	used for plants ar	nd bacteria.)	
4.	Water rinse 3 changes		37° C	60 sec. ea.	10 min. ea.
NOT	E: If using LR White	, Acetone or	PO can not	be used, only E	тон!
5.	Dehydration	50%	45° C	60 sec.	10 min.
	(Using either:	70%	45° C	60 sec.	10 min.
	ETOH, Acetone	80%	45° C	60 sec.	10 min.
	or Acetonitrile)	90%	45° C	60 sec.	10 min.
	2 changes	100%	45° C	60 sec. ea	10 min. ea
NOT	E: Separate the SEI	M for CPD or	HMDS from	the TEM samp	les at this stage if
nece	essary.				
(ETOH : Acetone OR		1:1	45° C	60 sec.	5 min.
ETOH : Acetonitrile OR 100% ETOH : PO)		100%	45° C	60 sec.	5 min.
6.	Infiltration — ETOH, Acetone, or PO:Resin				
	Pla	nt 3:1	50° C	15 min.	30 min.
		2:1	50° C	5 min.	1 hr.
		1:1	50° C	15 min.	1 hr.
	100% resin.	2 changes	50° C	15 min. ea.	1 hr. ea
7.	Embed in capsules and polymerize over night at 70° C				
8.	Or store in freezer in 100% resin until time for Embedment.				





Bone Marrow: Transmission electron microscope image of a thin section cut through an area of bone marrow area near the cartilage/bone interface in a mouse kneecap. Image shows small opening in the thin endothelium of the vascular sinus wall, where a blood cell is crossing the thin vascular sinus wall and into the sinus lumen.

Louisa Howard, Dartmouth College.

Buffers	
Phosphate	19340-72
Cacodylate	11652
Fixatives	
Karnovsky's Fixative EM Grade	15732-10
2% OsO ₄ Aqueous Solution	19172
Dehydrants	
Ethanol	15055
Acetone	10015
Propylene Oxide	20401
Acetonitrile	10020
Resins	
Embed 812	14120
Spurr's	14300
LR White	14383











