

Contact us today!

1560 Industry Road Hatfield, PA 19440
P.O. Box 550
TEL: 215-412-8400 FAX: 215-412-8450
TOLL FREE: 1-800-523-5874
EMAIL: sgkcck@aol.com
WEB: www.emsdiasum.com

INSTRUCTION MANUAL
CAT. 87023 Series
EMS MyBath™ Water Bath



About the EMS MyBath™ Water Bath

MyBath is a micro-processor controlled thermostatic water bath and is widely used for petroleum, biochemical and scientific research. The instrument is designed for accurate temperature control between ambient +5°C up to 99.9°C.

Product Specifications

Capacity	4L
Temperature range	Ambient +5°C to 99.9°C
Temperature accuracy	+/-0.2°C (37°C)
Temperature uniformity	+/-0.2°C (37°C)
Internal dimensions	w x d x h = 13.5 x 24 x 14 cm (5.3 x 9.4 x 5.5")
External dimensions	w x d x h = 19.5 x 36 x 24 cm (7.5 x 14 x 9.5")
Weight	6.5 kg (15 lbs.)
Electrical	100-120VAC or 220-240VAC, 50/60Hz, 8A/5A
Operating Environment	10°C to 30°C, 70%RH
Warranty	1 year

Product Set-Up

Place the water bath on a clean, level and stable surface. Make sure that the powerswitch is in the OFF position. Plug the power cord into the back of the unit, then plug the electrical cord into a properly grounded outlet of the correct voltage (shown on the serial number label). Prior to powering up the instrument, fill the tank with water (at least 50%).

WARNING: Before heating, ALWAYS ensure that the water tank is at least 50% full.

Product Operation

Once the tank is filled with liquid, power on the instrument (power switch is located on the back of the instrument). To set the desired temperature, press the "Temp." button and the current set temperature will flash on the display. Now use the up and down arrows to set the desired temperature. When the instrument is in the heating state, the hot indication light will illuminate.

WARNING: NEVER touch the inner chamber when the hot indication light is illuminated.

Recalibrating

Any adjustment to the temperature calibration can be done through the following:

1. Set the desired temperature to perform the calibration.
2. After this desired temperature has been reached, insert a thermometer and record the temperature.
3. Press the "Cal." key and the displayed temperature will flash.
4. Now use the up and down arrows to set the recorded temperature on the display. The calibration is stored after 8 seconds.

Electron Microscopy Sciences

1560 Industry Road Hatfield, PA 19440 P.O. Box 550
 TEL: 215-412-8400 FAX: 215-412-8450 TOLL FREE: 1-800-523-5874
 EMAIL: sgkcck@aol.com WEB: www.emsdiasum.com

Care and Maintenance

No routine maintenance is required other than to keep the unit clean. Cleaning can be done with a damp cloth. Avoid the use of solvents as they may attack the product housing.

Electron Microscopy Sciences

1560 Industry Road Hatfield, PA 19440 P.O. Box 550
TEL: 215-412-8400 FAX: 215-412-8450 TOLL FREE: 1-800-523-5874
EMAIL: sgkcck@aol.com WEB: www.emsdiasum.com