

INSTRUCTION MANUAL

CAT. 61040-12, 61040-20, 61040-EU, 61040-BR, 61040-SW, 61040-AU, 61040-10

Microplate Genie – Microplate Mixer – Analog



Electron Microscopy Sciences
1560 Industry Road Hatfield, PA 19440
TEL: 215-412-8400 FAX: 215-412-8450 TOLL FREE: 1-800-523-5874
EMAIL: sgkcck@aol.com WEB: www.emsdiasum.com

Caution:

CAUTION – Do not lift the Microplate Genie by the tray.

CAUTION – Failure to follow these operating instructions can compromise the user's safety.

CAUTION – Unplug from power before cleaning. Do not immerse.

CAUTION – Do not use Microplate Genie in hazardous atmospheres or with hazardous materials

CAUTION – Do not use Microplate Genie to mix flammable materials.

Care and handling:

Avoid wetting the Microplate Genie or exposing it to unnecessary fumes. After unplugging, you may wash the finish with water and soap or detergents with the use of a cloth or sponge. For best cleaning results, clean any blotting or spills immediately. For replacement parts, contact Electron Microscopy Sciences.

Electron Microscopy Sciences
1560 Industry Road Hatfield, PA 19440
TEL: 215-412-8400 FAX: 215-412-8450 TOLL FREE: 1-800-523-5874
EMAIL: sgkcck@aol.com WEB: www.emsdiasum.com