

► QUANTIFOIL® Holey Carbon Grids for cryo-TEM (continued)

Click on a product to add it to the cart

Ordering Information

Hole Shape	Grid Type	Hole Size	Period	Cat. # 200 Mesh	Cat. # 300 Mesh	Cat. # 400 Mesh	Qty.	
R 3.5/1	Copper	3.5µm	4.5µm	Q225CR-35	Q325CR-35	Q425CR-35	50/pk	
				Q250CR-35	Q350CR-35	Q450CR-35	100/pk	
	Nickel	3.5µm	4.5µm	Q225NR-35	Q325NR-35	Q425NR-35	50/pk	
				Q250NR-35	Q350NR-35	Q450NR-35	100/pk	
	Gold	3.5µm	4.5µm	Q225AR-35	Q325AR-35	Q425AR-35	50/pk	
				Q250AR-35	Q350AR-35	Q450AR-35	100/pk	
R 1/4	Copper	1µm	5µm	Q225CR-14	Q325CR-14	Q425CR-14	50/pk	
				Q250CR-14	Q350CR-14	Q450CR-14	100/pk	
	Nickel	1µm	5µm	Q225NR-14	Q325NR-14	Q425NR-14	50/pk	
				Q250NR-14	Q350NR-14	Q450NR-14	100/pk	
	Gold	1µm	5µm	Q225AR-14	Q325AR-14	Q425AR-14	50/pk	
				Q250AR-14	Q350AR-14	Q450AR-14	100/pk	
R 0.6/1	Copper	0.6µm	1.6µm	Q225CR-06	Q325CR-06	Q425CR-06	50/pk	
				Q250CR-06	Q350CR-06	Q450CR-06	100/pk	
	Nickel	0.6µm	1.6µm	Q225NR-06	Q325NR-06	Q425NR-06	50/pk	
				Q250NR-06	Q350NR-06	Q450NR-06	100/pk	
	Gold	0.6µm	1.6µm	Q225AR-06	Q325AR-06	Q425AR-06	50/pk	
				Q250AR-06	Q350AR-06	Q450AR-06	100/pk	
R 5/20	Copper	5µm	25µm	Q225CR-520	Q325CR-520	Q425CR-520	50/pk	
				Q250CR-520	Q350CR-520	Q450CR-520	100/pk	
	Nickel	5µm	25µm	Q225NR-520	Q325NR-520	Q425NR-520	50/pk	
				Q250NR-520	Q350NR-520	Q450NR-520	100/pk	
	Gold	5µm	25µm	Q225AR-520	Q325AR-520	Q425AR-520	50/pk	
				Q250AR-520	Q350AR-520	Q450AR-520	100/pk	
R 2/4	Copper	2µm	6µm	Q225CR-4	Q325CR4	Q425CR4	25/pk	
				Q250CR-4	Q350CR4	Q450CR4	50/pk	
				Q2100CR-4	Q3100CR4	Q4100CR4	100/pk	
	Nickel	2µm	6µm	Q225NR-4	Q325NR4	Q425NR4	25/pk	
				Q250NR-4	Q350NR4	Q450NR4	50/pk	
				Q2100NR-4	Q3100NR4	Q4100NR4	100/pk	
	Gold	2µm	6µm	Q225AR-4	Q325AR4	Q425AR4	25/pk	
				Q250AR-4	Q350AR4	Q450AR4	50/pk	
				Q2100AR-4	Q3100AR4	Q4100AR4	100/pk	
	R 2/2	Copper	2µm	4µm	Q225CR2	Q325CR2	Q425CR2	25/pk
					Q250CR2	Q350CR2	Q450CR2	50/pk
					Q2100CR2	Q3100CR2	Q4100CR2	100/pk
Nickel		2µm	4µm	Q225NR2	Q325NR2	Q425NR2	25/pk	
				Q250NR2	Q350NR2	Q450NR2	50/pk	
				Q2100NR2	Q3100NR2	Q4100NR2	100/pk	
Gold		2µm	4µm	Q225AR2	Q325AR2	Q425AR2	25/pk	
				Q250AR2	Q350AR2	Q450AR2	50/pk	
				Q2100AR2	Q3100AR2	Q4100AR2	100/pk	
R 2/1		Copper	2µm	3µm	Q225CR1	Q325CR1	Q425CR1	25/pk
					Q250CR1	Q350CR1	Q450CR1	50/pk
					Q2100CR1	Q3100CR1	Q4100CR1	100/pk
	Nickel	2µm	3µm	Q225NR1	Q325NR1	Q425NR1	25/pk	
				Q250NR1	Q350NR1	Q450NR1	50/pk	
				Q2100NR1	Q3100NR1	Q4100NR1	100/pk	
	Gold	2µm	3µm	Q225AR1	Q325AR1	Q425AR1	25/pk	
				Q250AR1	Q350AR1	Q450AR1	50/pk	
				Q2100AR1	Q3100AR1	Q4100AR1	100/pk	
	R 1.2/1.3	Copper	~1.2µm	~2.5µm	Q225CR1.3	Q325CR1.3	Q425CR1.3	25/pk
					Q250CR1.3	Q350CR1.3	Q450CR1.3	50/pk
					Q2100CR1.3	Q3100CR1.3	Q4100CR1.3	100/pk
Nickel		~1.2µm	~2.5µm	Q225NR1.3	Q325NR1.3	Q425NR1.3	25/pk	
				Q250NR1.3	Q350NR1.3	Q450NR1.3	50/pk	
				Q2100NR1.3	Q3100NR1.3	Q4100NR1.3	100/pk	
S 7/2	Copper	7 x 7µm	9µm	Q225CS7	Q325CS7	Q425CS7	25/pk	
				Q250CS7	Q350CS7	Q450CS7	50/pk	
				Q2100CS7	Q3100CS7	Q4100CS7	100/pk	
	Nickel	7 x 7µm	9µm	Q225NS7	Q325NS7	Q425NS7	25/pk	
				Q250NS7	Q350NS7	Q450NS7	50/pk	
				Q2100NS7	Q3100NS7	Q4100NS7	100/pk	
	Gold	7 x 7µm	9µm	Q225AS7	Q325AS7	Q425AS7	25/pk	
				Q250AS7	Q350AS7	Q450AS7	50/pk	
				Q2100AS7	Q3100AS7	Q4100AS7	100/pk	
	Multi A	Copper	—	—	Q225CMA	Q325CMA	Q425CMA	25/pk
					Q250CMA	Q350CMA	Q450CMA	50/pk
					Q2100CMA	Q3100CMA	Q4100CMA	100/pk
Nickel		—	—	Q225NMA	Q325NMA	Q425NMA	25/pk	
				Q250jNMA	Q350NMA	Q450NMA	50/pk	
				Q2100NMA	Q3100NMA	Q4100NMA	100/pk	
Hole Shape	Grid Type	Hole Size	Description	Cat. #	Pack			
HEX 15	Copper 2	6µm	0.5 x 2mm slot grids	Q225CR-HEX	50/pk			
				Q250CR-HEX	100/pk			