

III Acrylic 1mm Rat Heart Matrix

Acrylic heart matrix for the dissection of rat heart.

These heart matrices are precisely machined to ensure reproducible 1 mm sections of rat hearts.

SPECIFICATIONS

Slice	1 mm
Cavity Height (A)	19.9 mm
Cavity Height (B)	12.7 mm
Cavity Depth	9.7 mm

ORDERING INFORMATION

Cat. No.	Description	Qty.
69084-H	Acrylic 1 mm Rat Heart Matrix	each

