


## Spinal Cord Matrices


These Spinal Cord Matrices are precisely machined to ensure reproducible 0.5 mm or 1 mm sections of mouse or rat spinal cords.



69084-C



69085-C



69085-CS

## SPECIFICATIONS

Cat. No.	Cavity Size
69084-C	4 mm W x 96 mm L x 5 mm D
69085-C	4 mm W x 50 mm L x 5 mm D
69085-CS	2.5 mm W x 32 mm L x 5 mm D

## ORDERING INFORMATION

Cat. No.	Description	Qty.
69084-C	Spinal Cord Coronal, 0.5 mm Rat	each
69085-C	Spinal Cord Coronal, 0.5 mm Mouse	each
69085-CS	Spinal Cord Coronal & Sagittal, 1 mm Mouse	each