


III Stainless Steel Brain & Heart Matrices with 0.5mm Thick Slices

When your rodent brain, heart, or other similar organ needs to be sliced as thinly as possible, these new stainless steel matrices with a 0.5mm channel will achieve these slices – coronal or Sagittal – with ease.

FEATURES

- Made of high-quality stainless steel
- Can be repeatedly used, frozen, autoclaved, and cleaned
- Produces fine slices
- Accurate sectioning
- Easy cleaning and easy handling



SPECIFICATIONS

Cat. No.	Dimensions (mm)					Cavity Depth
	A	B	C	D	E	
69090-C	3.18	11.1	8.73	19.1	12.2	7.4
69090-S	3.18	11.1	8.73	19.1	12.2	7.4
69091-C	4.76	15.9	12.7	36.6	23.8	7.61
69091-S	4.76	15.9	12.7	36.6	23.8	7.61
69092-C	4.76	19.8	14.7	36.6	24.7	10.91
69092-S	4.76	19.5	14.7	36.6	24.7	10.91
69095-H	8.3	12.1	-	-	-	4.8
69096-H	12.7	19.8	-	-	-	9.6

ORDERING INFORMATION

Cat. No.	Description	Qty.
69090-C	Mouse Coronal (40-75g)	each
69090-S	Mouse Sagittal (40-75g)	each
69091-C	Small Rat Coronal (175-300g)	each
69091-S	Small Rat Sagittal (175-300g)	each
69092-C	Large Rat Coronal (300-600g)	each
69092-S	Large Rat Sagittal (300-600g)	each
69095-H	Mouse Heart	each
69096-H	Rat Heart	each