

III Compound Brightfield Microscopes

EMS keeps prices low, quality high. Our microscopes virtually eliminate down time and our optics deliver brilliant images. Your purchase will also be protected by our industry leading 5 year limited warranty.

Choose from the following options for the right equipment to meet your needs:

FEATURES

- Illumination: the popular rechargeable LED, or 20W halogen bulb or an incandescent bulb
- Monocular, Binocular or Teaching heads. Also, video systems are available to project images on a screen
- Three objective system, CP135-M3, has 4X, 10X and 40X
- Four objective system includes the 100X oil



SPECIFICATIONS

DIN Objective Magnification	Numerical Aperture (mm)	Field of View w/10X E.P. (mm)	Working Distance Brightfield DIN Achromatic (mm)
4X	0.1	4.5	27
10X	0.25	1.8	6.6
40X	0.65	0.45	0.57
100X Retractable (Oil)	1.25	0.18	0.19

ORDERING INFORMATION

Cat. No.	Description	Qty.
CP135-M3	Monocular, Education, 3 Achromat objectives	each
CP135-M4	Monocular, Education, 4 Achromat objectives	each
CP135-B4	Binocular, Education, 4 Achromat objectives	each
CP135-T4	Teaching	each